

# **Biotinylated Anti-Murine IL-5 Antibody**

Catalog # ABG10289

# **Specification**

## **Biotinylated Anti-Murine IL-5 Antibody - Product Information**

Application WB, E
Reactivity Mouse
Host Rabbit
Clonality Polyclonal

### Biotinylated Anti-Murine IL-5 Antibody - Additional Information

### **Preparation**

Produced from sera of rabbits immunized with highly pure recombinant Murine IL-5. Anti-Murine IL-5 specific antibody was purified by affinity chromatography and then biotinylated.

### WesternBlot

To detect Murine IL-5 by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 mg/ml. When used in conjunction with compatible development reagents, the detection limit for recombinant Murine IL-5 is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions. <span style="font-size: 12pt; font-family: &apos;times new roman&apos;, serif;"></span>

# Sandwich

To detect Murine IL-5 by sandwich ELISA (using 100 ml/well) a concentration of  $0.25 - 1.0 \mu g/ml$  of this antibody is required. This biotinylated polyclonal antibody, in conjunction with BioGems' Polyclonal Anti- Murine IL-5 (61-005P) as a capture antibody, allows the detection of at least  $0.2 - 0.4 \mu g/mell$  of recombinant Murine IL-5.<span style="font-size: 12pt; font-family: &apos;times new roman&apos;, serif;"></span>

### **Direct**

To detect Murine IL-5 by direct ELISA (using 100 ml/well) a concentration of approximately 1.0  $\mu$ g/ml of this antibody is required. This biotinylated polyclonal antibody allows the detection of at least 0.2 – 0.4 ng/well of recombinant Murine IL-5.

### **Formulation**

A sterile filtered antibody solution was lyophilized from PBS, pH 7.2.

### Reconstitution

Centrifuge vial prior to opening. Reconstitute in sterile PBS containing 0.1% BSA to a concentration of 0.1-1.0 mg/ml.

# **Storage**

-20°C

### **Precautions**

Biotinylated Anti-Murine IL-5 Antibody is for research use only and not for use in diagnostic or therapeutic procedures.





Tel: 858.875.1900 Fax: 858.875.1999

# **Biotinylated Anti-Murine IL-5 Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

**Biotinylated Anti-Murine IL-5 Antibody - Images**