

H3K27ac polyclonal antibody Rabbit Polyclonal Antibody Catalog # ABV11401

### **Specification**

## H3K27ac polyclonal antibody - Product Information

Application	
Primary Accession	
Host	
Clonality	
Isotype	
Calculated MW	

CHIP, DB, E, WB P68431 Rabbit Polyclonal Rabbit IgG 15404

#### H3K27ac polyclonal antibody - Additional Information

#### Gene ID 8350;8351;8352;8353;8354;8355;8356;8357;8358;8968

Positive Control

Western blot: Whole cell and histone extracts, IF: HeLa cells, ELISA: Antigen, ChIP: HeLa cells, Dot blot: Histone Peptides. WB: 1:1000, Dot Blot: 1:25,000, ChIP: 1-2 µl/ChIP, IF: 1:500, ELISA: 1:100.

Application & Usage

Other Names Histone 3

Target/Specificity H3K27ac

Antibody Form Liquid

Appearance Colorless liquid

Formulation In PBS with 0.05% sodium azide and 0.05% ProClin 300.

Handling The antibody solution should be gently mixed before use.

Reconstitution & Storage -20 °C

**Background Descriptions** 

### Precautions

H3K27ac polyclonal antibody is for research use only and not for use in diagnostic or therapeutic procedures.



# H3K27ac polyclonal antibody - Protein Information

Name H3C1 (<u>HGNC:4766</u>)

Synonyms H3FA, HIST1H3A

Function

Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling.


Cellular Location Nucleus. Chromosome.

## H3K27ac polyclonal antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.


- Western Blot
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- <u>Cell Culture</u>

# H3K27ac polyclonal antibody - Images




ChIP assays were performed using HeLa cells and the antibody and optimized PCR primer sets for qPCR. A titration of the antibody consisting of 1, 2, 5, and 10  $\mu$ I per ChIP experiment was analysed. IgG (2  $\mu$ g/IP) was used as negative control. The Fig shows the recovery, expressed as a % of input (the relative amount of IP DNA compared to input DNA after qPCR analysis).






A Dot Blot analysis was performed to test the cross reactivity of the antibody with peptides containing other modifications of histone H4 and H3. 100 to 0.2 pmol of the peptide containing the respective histone modification were spotted on a membrane. The Fig shows a high specificity of the antibody for the modification of interest.



To determine the titer, an ELISA was performed using a serial dilution of the antibody. The antigen used was a peptide containing the histone modification of interest. By plotting the absorbance against the antibody dilution the titer of the antibody was estimated to be 1:4,000.



Western blot was performed on whole cell (40  $\mu$ g, lane 1) and histone extracts (15  $\mu$ g, lane 2) from HeLa cells, and on 1  $\mu$ g of recombinant histone H2A, H2B, H3 and H4 (lane 3, 4, 5 and 6, respectively) using the antibody. The position of the protein of interest is shown on the right, the marker (in kDa) is shown on the left.

# H3K27ac polyclonal antibody - Background

Histones are the main constituents of the protein part of chromosomes of eukaryotic cells. They are rich in the amino acids arginine and lysine and have been greatly conserved during evolution.



Histones pack the DNA into tight masses of chromatin. Two core histones of each class H2A, H2B, H3 and H4 assemble and are wrapped by 146 base pairs of DNA to form one octameric nucleosome. Histone tails undergo numerous post-translational modifications, which either directly or indirectly alter chromatin structure to facilitate transcriptional activation or repression or other nuclear processes. In addition to the genetic code, combinations of the different histone modifications reveal the so-called "histone code". Histone methylation and demethylation is dynamically regulated by respectively histone methyl transferases and histone demethylases.