

CD74 (B-Cell Marker) Antibody - With BSA and Azide

Mouse Monoclonal Antibody [Clone CLIP/1133] **Catalog # AH12809**

Specification

CD74 (B-Cell Marker) Antibody - With BSA and Azide - Product Information

Application ,1,3,4, **Primary Accession** P044233 Other Accession 972, 436568 Reactivity Human Host Mouse Clonality **Monoclonal** Isotype

Mouse / IgG1, kappa

Calculated MW 33-41kDa KDa

CD74 (B-Cell Marker) Antibody - With BSA and Azide - Additional Information

Storage

Store at 2 to 8°C.Antibody is stable for 24 months.

Precautions

CD74 (B-Cell Marker) Antibody - With BSA and Azide is for research use only and not for use in diagnostic or therapeutic procedures.

CD74 (B-Cell Marker) Antibody - With BSA and Azide - Protein Information


CD74 (B-Cell Marker) Antibody - With BSA and Azide - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

CD74 (B-Cell Marker) Antibody - With BSA and Azide - Images





CD74 (B-Cell Marker) Antibody - With BSA and Azide - Background

Recognizes proteins of 33, 35 and 41kDa, which are identified as various isoforms of CD74. Its epitope is localized in the extracellular domain of CD74. CD74 is a type II transmembrane protein which binds to the peptide binding groove of newly synthesized MHC class II alpha/beta heterodimers and prevents their premature association with endogenous polypeptides. The CD74 molecule plays a critical role in the presentation of peptides, by the MHC class II antigens, to CD4 positive lymphocytes. CD74 is expressed on MHC class II positive cells including B cells, a subset of activated T cells, monocytes, and dendritic cells and by various types of carcinomas. CD74 is expressed primarily by antigen presenting cells, such as B-lymphocytes (from before the pre-B cell stage to before the plasma cell stage), macrophages, and monocytes, and many epithelial cells. Anti-CD74 stains predominantly germinal center lymphocytes and B-cell lymphomas, but rarely T-cell lymphomas. Anti-CD74 has been shown to be useful in differentiating atypical fibroxanthoma (-) from malignant fibrous histiocytoma (+).

CD74 (B-Cell Marker) Antibody - With BSA and Azide - References

Wraight, C.J. et al. (1990) Human major histocompatibility complex class II invariant chain is expressed on the cell surface. J Biol Chem. 265:5787