

Swine IgG Peroxidase

Catalog # ASR1082

Specification

Swine IgG Peroxidase - Product Information

Description **SWINE IgG whole molecule Peroxidase**

conjugated

Conjugate Peroxidase (Horseradish) **Physical State**

Lyophilized

Host Isotype IgG Buffer

0.02 M Potassium Phosphate, 0.15 M

Sodium Chloride, pH 7.2

Species of Origin **Swine** Reconstitution Volume 1.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Gentamicin Sulfate. Do NOT

add Sodium Azide!

Swine IgG Peroxidase - Additional Information

Shipping Condition

Ambient

Purity

This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase, anti-Swine IgG and anti-Swine Serum.

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Swine IgG Peroxidase - Protein Information

Swine IgG Peroxidase - Protocols

Provided below are standard protocols that you may find useful for product applications.





Tel: 858.875.1900 Fax: 858.875.1999

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- <u>Immunoprecipitation</u>
- Flow Cytomety
- Cell Culture

Swine IgG Peroxidase - Images