

Dog IgG Fab Rhodamine

Catalog # ASR1264

Specification

Dog IgG Fab Rhodamine - Product Information

Description DOG IgG F(ab) fragment Rhodamine

conjugated

Conjugate Rhodamine (TRITC)

FP Value 1.2 moles Rhodamine (TRITC) per mole of

Dog IgG F(ab)
Lyophilized

Physical State
Host Isotype
Buffer
Lyophilized
IgG F(ab)
0.02 M Pot

0.02 M Potassium Phosphate, 0.15 M

Sodium Chloride, pH 7.2

Dog 1.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Sodium Azide

Dog IgG Fab Rhodamine - Additional Information

Shipping Condition

Species of Origin

Reconstitution Volume

Ambient

Purity

This product was prepared from normal serum by delipidation, salt fractionation, ion exchange chromatography followed by pepsin digestion and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Dog IgG, anti-Dog IgG F(ab')2 and anti-Dog Serum. No reaction was observed against anti-Dog IgG F(c) or anti-Pepsin.

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Dog IgG Fab Rhodamine - Protein Information

Dog IgG Fab Rhodamine - Protocols





Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- <u>Immunofluorescence</u>
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Dog IgG Fab Rhodamine - Images

Dog IgG Fab Rhodamine - Background

This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.