

Horse IgG F(ab')2 Fluorescein

Catalog # ASR1277

Specification

Horse IgG F(ab')2 Fluorescein - Product Information

Description HORSE IgG F(ab')2 fragment

Conjugate Fluorescein (FITC)

FP Value 2.3 moles Fluorescein (FITC) per mole of

Horse IgG F(ab')2

Physical State
Host Isotype
Lyophilized
IgG F(ab')2

Buffer 0.01 M Sodium Phosphate, 0.15 M Sodium

Chloride, pH 7.2

Species of Origin
Reconstitution Volume
1.0 mL

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer 10 mg/mL Bovine Serum Albumin (BSA) -

Immunoglobulin and Protease free

Preservative 0.01% (w/v) Thimerosal

Horse IgG F(ab')2 Fluorescein - Additional Information

Shipping Condition

Ambient

Purity

This product was prepared from normal serum by delipidation, salt fractionation and ion change chromatography followed by pepsin digestion and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Fluorescein, anti-Horse IgG, anti-Horse IgG F(ab')2 and anti-Horse Serum. No reaction was observed against anti-Horse IgG F(c) or anti-Pepsin.

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Horse IgG F(ab')2 Fluorescein - Protein Information

Horse IgG F(ab')2 Fluorescein - Protocols





Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Horse IgG F(ab')2 Fluorescein - Images

Horse IgG F(ab')2 Fluorescein - Background

This product is designed for immunofluorescence microscopy, fluorescence based plate assays (FLISA) and fluorescent western blotting. This product is also suitable for multiplex analysis, including multicolor imaging, utilizing various commercial platforms.