

Anti-Guinea Pig IgG (H&L) Pre-Adsorbed Secondary Antibody Goat Polyclonal, Unconjugated Catalog # ASR1391

Specification

Anti-Guinea Pig IgG (H&L) Pre-Adsorbed Secondary Antibody - Product Information

Description

Host Conjugate Target Species Clonality Application Application Note

Physical State Host Isotype Target Isotype Buffer

Immunogen Stabilizer Preservative Anti-GUINEA PIG IgG (H&L) (GOAT) Antibody (Min X By Ch Gt Ham Hs Hu Ms Rb **Rt & Sh Serum Proteins)** Goat Unconjugated **Guinea** Pig Polyclonal ,1,10,15, ELISA 1:50,000;Western Blot 1:2,000-1:10,000;Immunochemistry 1:1.000-1:5.000 Liquid (sterile filtered) laG IgG (H&L) 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 Guinea Pig IgG whole molecule None 0.01% (w/v) Sodium Azide

Anti-Guinea Pig IgG (H&L) Pre-Adsorbed Secondary Antibody - Additional Information

Shipping Condition Wet Ice

Purity

This product was prepared from monospecific antiserum by immunoaffinity chromatography using Guinea Pig IgG coupled to agarose beads followed by solid phase adsorption(s) to remove any unwanted reactivities. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Goat Serum, Guinea Pig IgG and Guinea Pig Serum. No reaction was observed against Bovine, Chicken, Goat, Hamster, Human, Mouse, Rabbit, Rat and Sheep Serum Proteins. Minimal cross reactivity was observed against the species listed above.

Storage Condition

Store vial at 4° C prior to opening. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.

Anti-Guinea Pig IgG (H&L) Pre-Adsorbed Secondary Antibody - Protein Information



Anti-Guinea Pig IgG (H&L) Pre-Adsorbed Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- <u>Western Blot</u>
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- <u>Cell Culture</u>

Anti-Guinea Pig IgG (H&L) Pre-Adsorbed Secondary Antibody - Images