

F(ab')2 Anti-Guinea Pig IgG F(ab')2 Secondary Antibody

Rabbit Polyclonal, Unconjugated Catalog # ASR3449

Specification

F(ab')2 Anti-Guinea Pig IgG F(ab')2 Secondary Antibody - Product Information

F(ab')2 Anti-GUINEA PIG IgG Description F(ab')2 (RABBIT) Antibody

Host **Rabbit** Unconjugated Conjugate **Target Species Guinea Pig** Clonality **Polyclonal**

Application ,1,2,10,

Application Note ELISA 1:20,000-1:100,000; Western Blot 1:2,000-1:10,000;Immunohistochemistry

1:1,000-1:5,000

Lyophilized **Physical State** Host Isotype IgG F(ab')2 IgG F(ab')2 Target Isotype

Buffer 0.02 M Potassium Phosphate, 0.15 M

Sodium Chloride, pH 7.2

Immunogen Guinea Pig IgG F(ab')2 fragment

Reconstitution Volume

Reconstitution Buffer Restore with deionized water (or

equivalent)

Stabilizer None

F(ab')2 Anti-Guinea Pig IgG F(ab')2 Secondary Antibody - Additional Information

Shipping Condition

Ambient

Purity

This product is a F(ab')2 fragment of IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation, ion exchange chromatography and pepsin digestion followed by chromatographic separation and extensive dialysis against the buffer stated above. Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum, Guinea Pig IgG, Guinea Pig IgG F(ab')2 and Guinea Pig Serum. No reaction was observed against Guinea Pig IgG F(c), anti-Rabbit IgG F(c) or anti-Pepsin.

Storage Condition

Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.

Precautions Note

This product is for research use only and is not intended for therapeutic or diagnostic applications.



F(ab')2 Anti-Guinea Pig IgG F(ab')2 Secondary Antibody - Protein Information

F(ab')2 Anti-Guinea Pig IgG F(ab')2 Secondary Antibody - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- <u>Immunohistochemistry</u>
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

F(ab')2 Anti-Guinea Pig IgG F(ab')2 Secondary Antibody - Images

F(ab')2 Anti-Guinea Pig IgG F(ab')2 Secondary Antibody - Background

F(ab')2 Antibody was generated by enzymatic cleavage and subsequent separation from the Fc fragment. Because of their smaller size, F(ab)2 fragments offer several advantages over intact antibodies for use in certain immunochemical techniques and experimental applications. F(ab)2 fragments penetrate into tissue samples and show better antigen recognition and signal generation in IHC. F(ab)2 fragments lack the Fc region and therefore do not bind Fc receptors which effectively lowers background staining. F(ab')2 Antibody is ideal for investigators who routinely perform flow cytometry, immunohistochemistry or IHC and other immunoassays.