

## **HIS Tag Antibody**

Mouse Monoclonal Antibody (Mab)
Catalog # AW5672

# **Specification**

## **HIS Tag Antibody - Product Information**

Application Host Clonality Calculated MW Isotype Antigen Source WB,E Mouse Monoclonal 45-50KD KDa Mouse IgG1 HUMAN

## **HIS Tag Antibody - Additional Information**

## **Antigen Region**

NA

#### **Dilution**

WB~~1:1000

#### **Target/Specificity**

Purified recombinant HIS-tagged fusion protein and poly-HIS peptide were used to produced this monoclonal antibody.

## Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

#### **Precautions**

HIS Tag Antibody is for research use only and not for use in diagnostic or therapeutic procedures.

#### **HIS Tag Antibody - Protein Information**


# **HIS Tag Antibody - Protocols**

Provided below are standard protocols that you may find useful for product applications.


- Western Blot
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture



## **HIS Tag Antibody - Images**



Anti-HIS Tag at 1:8000 dilution + 12-tag lysate Lysates/proteins at 20ng per lane. Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size : 48 kDa Blocking/Dilution buffer: 5% NFDM/TBST.




All lanes: Anti-HIS Tag at 1:1000 dilution (upper) or GAPDH (lower) Lane 1: 293T/17 transfected with 12tag lysate (10ug) Lane 2: Non-transfected 293T/17 lysate (10ug) Secondary Goat Anti-mouse IgG, (H+L), Peroxidase conjugated at 1/10000 dilution. Predicted band size: 45-50 kDa Blocking/Dilution buffer: 5% NFDM/TBST.

## **HIS Tag Antibody - Background**

Epitope tags consisting of short sequences recognized by well-characterized monoclonal antibodies have been widely used in the study of protein expression in various systems. The 6xHIS tag (HHHHHH), recognized by the monoclonal antibody clone 6AT18 provides an established example of this application. 6xHIS-tagged fusion proteins are easily purified from cell lysates by affinity chromatography using Nickel-Sepharose resin. Abgent's anti-6xHIS monoclonal antibody provides a simple solution to detect the expression of HIS-tagged fusion proteins in cells.

#### **HIS Tag Antibody - References**





Hochuli E, Doebeli H, and Schacher A. New metal chelate absorbent selective for proteins and peptides containing neighboring histidine residues. J. Chromatogr. 1987;411:177-184.