

USP22 Antibody (C-term) Blocking Peptide
Synthetic peptide
Catalog # BP2148b**Specification**

USP22 Antibody (C-term) Blocking Peptide - Product InformationPrimary Accession [Q9UPT9](#)**USP22 Antibody (C-term) Blocking Peptide - Additional Information****Gene ID** 23326**Other Names**

Ubiquitin carboxyl-terminal hydrolase 22, Deubiquitinating enzyme 22, Ubiquitin thioesterase 22, Ubiquitin-specific-processing protease 22, USP22, KIAA1063, USP3L

Target/Specificity

The synthetic peptide sequence used to generate the antibody [AP2148b](/product/products/AP2148b) was selected from the C-term region of human USP22. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

Format

Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

Precautions

This product is for research use only. Not for use in diagnostic or therapeutic procedures.

USP22 Antibody (C-term) Blocking Peptide - Protein Information**Name** USP22**Synonyms** KIAA1063, USP3L**Function**

Histone deubiquitinating component of the transcription regulatory histone acetylation (HAT) complex SAGA. Catalyzes the deubiquitination of both histones H2A and H2B, thereby acting as a coactivator. Recruited to specific gene promoters by activators such as MYC, where it is required for transcription. Required for nuclear receptor-mediated transactivation and cell cycle progression.

Cellular Location

Nucleus {ECO:0000250|UniProtKB:Q5DU02}.

Tissue Location

Moderately expressed in various tissues including heart and skeletal muscle, and weakly expressed in lung and liver

USP22 Antibody (C-term) Blocking Peptide - Protocols

Provided below are standard protocols that you may find useful for product applications.

- [Blocking Peptides](#)

USP22 Antibody (C-term) Blocking Peptide - Images**USP22 Antibody (C-term) Blocking Peptide - Background**

Modification of target proteins by ubiquitin participates in a wide array of biological functions. Proteins destined for degradation or processing via the 26 S proteasome are coupled to multiple copies of ubiquitin. However, attachment of ubiquitin or ubiquitin-related molecules may also result in changes in subcellular distribution or modification of protein activity. An additional level of ubiquitin regulation, deubiquitination, is catalyzed by proteases called deubiquitinating enzymes, which fall into four distinct families. Ubiquitin C-terminal hydrolases, ubiquitin-specific processing proteases (USPs),¹ OTU-domain ubiquitin-aldehyde-binding proteins, and Jab1/Pad1/MPN-domain-containing metallo-enzymes. Among these four families, USPs represent the most widespread and represented deubiquitinating enzymes across evolution. USPs tend to release ubiquitin from a conjugated protein. They display similar catalytic domains containing conserved Cys and His boxes but divergent N-terminal and occasionally C-terminal extensions, which are thought to function in substrate recognition, subcellular localization, and protein-protein interactions.

USP22 Antibody (C-term) Blocking Peptide - References

Kikuno, R., et al., DNA Res. 6(3):197-205 (1999).