

Thrombin, Active, Porcine blood recombinant protein

Activated Factor IIa Catalog # PBV11197r

Specification

Thrombin, Active, Porcine blood recombinant protein - Product info

Primary Accession <u>F1SIB1</u>

Calculated MW 33.5 kDa KDa

Thrombin, Active, Porcine blood recombinant protein - Additional Info

Gene Symbol F2

Other Names Activated Factor IIa

Gene Source Porcine

Source Porcine blood Assay&Purity SDS-PAGE; ≥98%

Assay2&Purity2 N/A;
Recombinant No

Results 20-1000 IU/mg of lyophilized powder

Target/Specificity

Thrombin

Application Notes

Reconstitute in sterile water (100 U/ml) with 0.9% NaCl. It forms a clear solution.

Format Lyophilized

Storage

-20°C; Sterile filtered and lyophilized without any salts.

Thrombin, Active, Porcine blood recombinant protein - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Thrombin, Active, Porcine blood recombinant protein - Images

Thrombin, Active, Porcine blood recombinant protein - Background





Tel: 858.875.1900 Fax: 858.875.1999

Thrombin enzyme (Activated Factor IIa) is an important clotting promoter that controls the transformation of soluble fibringgen to insoluble active fibrin strands. Thrombin is a coagulation protein and a serine protease (EC 3.4.21.5) that catalyzes many coagulation-related reactions. Thrombin triggers factor-XI, factor-V, Factor-XIII and factor-VIII. Thrombin endorses platelet activation, using activation of protease-activated receptors on the platelet. As a result of its high proteolytic specificity, thrombin has become an important biochemical protein. The thrombin cleavage site (Leu-Val-Pro-Arg-Gly-Ser) is widely used in linker regions of recombinant fusion protein constructs. After the purification of the fusion protein, thrombin is used to cleave between the Arginine and Glycine residues of the cleavage site, efficiently removing the purification tag from the protein of interest with a high degree of specificity.