

Thrombin, Active, Bovine Plasma (High Activity) recombinant protein

Activated Factor IIa Catalog # PBV11199r

Specification

Thrombin, Active, Bovine Plasma (High Activity) recombinant protein - Product info

Primary Accession P00735
Calculated MW 37 kDa KDa

Thrombin, Active, Bovine Plasma (High Activity) recombinant protein - Additional Info

Gene ID 280685
Gene Symbol F2

Other Names
Activated Factor IIa

Gene Source Bovine

Source Bovine plasma
Assay&Purity SDS-PAGE; ≥98%

Assay2&Purity2 N/A;
Recombinant No

Results >1500 units/mg powder.

Target/Specificity

Thrombin

Application Notes

Reconstitute in sterile water (100 U/ml) with 0.9% NaCl. It forms a clear solution.

Format Lyophilized

, ,

Storage

-20°C; Sterile filtered and lyophilized with Mannitol and Sodium Chloride.

Thrombin, Active, Bovine Plasma (High Activity) recombinant protein - Protocols

Provided below are standard protocols that you may find useful for product applications.

- Western Blot
- Blocking Peptides
- Dot Blot
- Immunohistochemistry
- Immunofluorescence
- Immunoprecipitation
- Flow Cytomety
- Cell Culture

Thrombin, Active, Bovine Plasma (High Activity) recombinant protein - Images



Thrombin, Active, Bovine Plasma (High Activity) recombinant protein - Background

Thrombin enzyme (Activated Factor IIa) is an important clotting promoter that controls the transformation of soluble fibrinogen to insoluble active fibrin strands. Thrombin is a coagulation protein and a serine protease (EC 3.4.21.5) that catalyzes many coagulation-related reactions. Thrombin triggers factor-XI, factor-V, Factor-XIII and factor-VIII. Thrombin endorses platelet activation, using activation of protease-activated receptors on the platelet. As a result of its high proteolytic specificity, thrombin has become an important biochemical protein. The thrombin cleavage site (Leu-Val-Pro-Arg-Gly-Ser) is widely used in linker regions of recombinant fusion protein constructs. After the purification of the fusion protein, thrombin is used to cleave between the Arginine and Glycine residues of the cleavage site, efficiently removing the purification tag from the protein of interest with a high degree of specificity.

Thrombin, Active, Bovine Plasma (High Activity) recombinant protein - References

Irwin D.M.,et al.J. Mol. Biol. 200:31-45(1988). McGillivray R.T.A.,et al.Biochemistry 23:1626-1634(1984). Magnusson S.,et al.(In) Hemker H.C., Veltkamp J.J. (eds.); Irwin D.M.,et al.Biochemistry 24:6854-6861(1985). Park C.H.,et al.Biochemistry 25:3977-3982(1986).