

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic
Synthetic Peptide
Catalog # SP3696a**

Specification

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Product
Information**

Sequence **Dnp-RPLALWRS-COOH**

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Additional
Information**

Format

Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

Precautions

This product is for research use only. Not for use in diagnostic or therapeutic procedures.

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Protein
Information**

MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Images