

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic  
Synthetic Peptide  
Catalog # SP3696b**

**Specification**

---

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Product  
Information**

Sequence **Dnp-RPLALWRS-COOH**

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Additional  
Information**

**Format**

Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

**Storage**

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C.

**Precautions**

This product is for research use only. Not for use in diagnostic or therapeutic procedures.

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Protein  
Information**

**MMP-7 Substrate (Matrilysin Substrate, PUMP-1 Substrate), Fluorogenic - Images**